

Nanotrap[®] Microbiome Combined; 10 mL Automated Protocol using NucleoMag[®] Kit and the KingFisher[™] Flex

Objective: This protocol uses Nanotrap Microbiome A Particles and Nanotrap Microbiome B Particles and Nanotrap Enhancement Reagent 3 to capture and concentrate microbes in environmental water samples. It is optimized for microbe capture from 10 mL samples and is compatible with MACHEREY-NAGEL NucleoMag DNA/RNA Water Kit. The automated script can process up to 24 samples at once and can be amended for the throughput in your lab.

Materials and equipment:

Sample Type	
Environmental Water Samples	
Concentration Reagent	Vendor
Nanotrap Microbiome A Particles	Ceres Nanosciences; SKU# 44202
Nanotrap Microbiome B Particles	Ceres Nanosciences; SKU# 65202
Nanotrap Enhancement Reagent 3 (ER3) ¹	Ceres Nanosciences; SKU# 10113
Extraction Kit	Vendor
NucleoMag DNA/RNA Water Extraction Kit	MACHEREY-NAGEL; REF 744220.1
Materials/Equipment	Vendor
KingFisher [™] Flex Purification System, KingFisher with 96 Deep-well Head	Thermo Fisher Scientific [™] ; Cat# 5400630
KingFisher [™] Flex 24 Deep Well head	Thermo Fisher Scientific [™] ; Cat# 24074440
KingFisher [™] Flex 24 Deep Well heating block	Thermo Fisher Scientific [™] ; Cat# 24075440
KingFisher [™] Flex 96 heating block	Thermo Fisher Scientific [™] ; Cat# 24075420
KingFisher [™] 24 deep-well plate (for Duo Prime, Flex and Presto)	Thermo Fisher Scientific [™] ; Cat# 95040470
KingFisher [™] 24 deep-well tip comb and plate (for Flex and Presto)	Thermo Fisher Scientific [™] ; Cat# 97002610
KingFisher [™] 96 deep-well plate, v-bottom, polypropylene (for Duo Prime, Flex and Presto)	Thermo Fisher Scientific [™] ; Cat# 95040450
KingFisher [™] 96 tip comb for deep-well magnets, 10 x 10 pcs/box (for Flex and Presto)	Thermo Fisher Scientific [™] ; Cat# 97002534
General Reagents	Vendor
Molecular grade water	VWR; 45001-044

¹ Precipitate can form in ER3 if stored below room temperature. Allow ER3 to return to room temperature to dissolve the precipitate (can invert 2-3 times to help resuspend, do not heat).

Capture and Extract Microbes using Nanotrap Microbiome Particles

Procedure:

1. Nanotrap Microbiome Combined NucleoMag KingFisher Flex Procedure-Part 1

1. *Prepare* “Sample Plate 1” and “Sample Plate 2”
 1. Invert environmental water sample 5 times to mix. After inverting, place on a flat surface for 45 seconds.
 2. Add 4,800 µL of environmental water sample to one well (one well per sample) of a new KingFisher 24 Well Deep Well Plate.
 3. Add another 4,800 µL of environmental water sample to the same well location on a second KingFisher™ 24 Well Deep Well Plate.
 - a) For example, if you loaded a sample into well A1 of the first plate, load the second volume of that sample into well A1 of the second plate.
 4. Add 50 µL of Nanotrap Enhancement Reagent 3 (ER3) Solution to each sample on the two KingFisher 24 Well Deep Well sample plates (100 µL total).
 5. Add 75 µL of Nanotrap® Microbiome A Particles to each sample on the two KingFisher 24 Well Deep Well sample plates (150 µL total).
 6. Add 75 µL of Nanotrap Microbiome B Particles to each sample on the two KingFisher 24 Well Deep Well sample plates (150 µL total).
2. *Prepare* “Lysis Plate”
 1. Add 500 µL of Lysis Buffer MWA1 to a new (the third) KingFisher 24 Well Deep Well Plate matching the number and location of the “Sample Plate” wells.
3. *Prepare* “Tip Plate”
 1. Insert a new tip comb into a new KingFisher 24 Well Deep Well Plate.
4. *Run NT Script (Request file at sales@ceresnano.com)*
 1. Run
NT_Microbiome_A_B_NucleoMag®_24_w_heat_Flex_10mL.bdz
 2. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
5. Once the protocol is completed, the “Lysis Plate” will contain lysate that is ready to proceed to Part 2 (***caution* sample may be hot**).

2. Nanotrap Microbiome Combined NucleoMag KingFisher Flex Procedure-Part 2

1. *Prepare* “NM Binding” Plate
 1. To a new KingFisher 96 Deep Well Plate, add 450 µL of the cleared lysate (NT lysate) from each well of the lysis plate used in “Part 1 step 5” of the

protocol. Keep track of which well contains which sample in this new bead binding plate.

2. Add 475 μ L of Binding Buffer MWA2 to each well in which lysate was added.
 3. Vortex the NucleoMag B-beads thoroughly and add 25 μ L to each well.
 - a) Note: Binding mix (MWA2 + B-beads) can be pre-mixed before their addition to the plate.
 2. Prepare "1st Wash MWA3" Plate
 1. Add 850 μ L of Wash Buffer MWA3 to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- "NM Binding" Plate wells.
 3. Prepare "2nd Wash MWA3" Plate
 1. Add 850 μ L of Wash Buffer MWA3 to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- "NM Binding" Plate wells.
 4. Prepare "3rd Wash MWA4" Plate
 1. Add 850 μ L of Wash Buffer MWA4 to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- "NM Binding" Plate wells.
 5. Prepare "Elution" Plate
 1. Add 100 μ L of Rnase-free water to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- "NM Binding" Plate wells.
 6. Prepare "Tip Plate"
 1. Insert the KingFisher 96 Deep Well Comb into a new KingFisher 96 Deep Well Plate
 7. *Run Extraction Kit Script (Request file at sales@ceresnano.com)*
 1. Run
NucleoMag_DNA_RNA_Water_CeresNanoTrap_Flex_Rev02.bdz
 2. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
3. Once the protocol is completed, the KingFisher 96-Elution Plate contains eluates that are ready for downstream analysis or can be stored at -80°C.
Note: Multiple freeze-thaw cycles may cause degradation.

Attachments: 2

KingFisher™ Flex

1. *NT_Microbiome_A_and_B_NucleoMag®_24_w_heat_Flex_10mL.bdz*
2. *NucleoMag_DNA_RNA_Water_CeresNanoTrap_Flex_Rev02.bdz*