

## Nanotrap<sup>®</sup> Microbiome Combined; 10 mL Automated Protocol using NucleoMag<sup>®</sup> Kit and the KingFisher<sup>™</sup> Apex

**Objective:** This protocol uses Nanotrap Microbiome A Particles and Nanotrap Microbiome B Particles and Nanotrap Enhancement Reagent 3 to capture and concentrate microbes in environmental water samples. It is optimized for microbe capture from 10 mL samples and is compatible with MACHEREY-NAGEL NucleoMag DNA/RNA Water Kit. The automated script can process up to 24 samples at once and can be amended for the throughput in your lab.

### Materials and equipment:

Sample Type	
Environmental Water Samples	
Concentration Reagent	Vendor
Nanotrap Microbiome A Particles	Ceres Nanosciences; SKU# 44202
Nanotrap Microbiome B Particles	Ceres Nanosciences; SKU# 65202
Nanotrap Enhancement Reagent 3 (ER3) <sup>1</sup>	Ceres Nanosciences; SKU# 10113
Extraction Kit	Vendor
NucleoMag DNA/RNA Water Extraction Kit	MACHEREY-NAGEL; REF 744220.1
Materials/Equipment	Vendor
KingFisher <sup>™</sup> Apex with 96 DW Head	Thermo Fisher Scientific <sup>™</sup> ; Cat# 5400930
KingFisher Apex 24 Combi head	Thermo Fisher Scientific; Cat# 24079940
KF Apex 96 KF heating block	Thermo Fisher Scientific; Cat# 24075920
KF Apex 24 DW heating block	Thermo Fisher Scientific; Cat# 24075940
KingFisher 24 Deep-well Plate, Barcoded	Thermo Fisher Scientific; Cat#95040470B
KingFisher 24 Deep-Well Tip Comb & Plate, Barcoded	Thermo Fisher Scientific; Cat#97002610B
KingFisher 96 Deep-well Plate, Barcoded	Thermo Fisher Scientific; Cat# 95040450B
KingFisher 96 Plate (200 µL), Barcoded	Thermo Fisher Scientific; Cat# 97002540B
KingFisher 96 Deep-well Tip Comb, Barcoded	Thermo Fisher Scientific; Cat# 97002534B
General Reagents	Vendor
Molecular grade water	VWR; 45001-044

<sup>1</sup> Precipitate can form in ER3 if stored below room temperature. Allow ER3 to return to room temperature to dissolve the precipitate (can invert 2-3 times to help resuspend, do not heat).

## Capture and Extract Microbes using Nanotrap Microbiome Particles

### Procedure:

#### 1. Nanotrap Microbiome Combined NucleoMag KingFisher Apex Procedure-Part 1

1. *Prepare* “Sample Plate 1” and “Sample Plate 2”
  1. Invert environmental water sample 5 times to mix. After inverting, place on a flat surface for 45 seconds.
  2. Add 4,800  $\mu\text{L}$  of environmental water sample to one well (one well per sample) of a new KingFisher 24 Well Deep Well Plate.
  3. Add another 4,800  $\mu\text{L}$  of environmental water sample to the same well location on a second KingFisher™ 24 Well Deep Well Plate.
    - a) For example, if you loaded a sample into well A1 of the first plate, load the second volume of that sample into well A1 of the second plate.
  4. Add 50  $\mu\text{L}$  of Nanotrap Enhancement Reagent 3 (ER3) Solution to each sample on the two KingFisher 24 Well Deep Well sample plates (100  $\mu\text{L}$  total).
  5. Add 75  $\mu\text{L}$  of Nanotrap® Microbiome A Particles to each sample on the two KingFisher 24 Well Deep Well sample plates (150  $\mu\text{L}$  total).
  6. Add 75  $\mu\text{L}$  of Nanotrap Microbiome B Particles to each sample on the two KingFisher 24 Well Deep Well sample plates (150  $\mu\text{L}$  total).
2. *Prepare* “Lysis Plate”
  1. Add 500  $\mu\text{L}$  of Lysis Buffer MWA1 to a new (the third) KingFisher 24 Well Deep Well Plate matching the number and location of the “Sample Plate” wells.
3. *Prepare* “Tip Plate”
  1. Insert a new tip comb into a new KingFisher 24 Well Deep Well Plate.
4. *Run NT Script* (Request file at [sales@ceresnano.com](mailto:sales@ceresnano.com))
  1. Run **NT\_Microbiome\_A\_and\_B\_NucleoMag®\_24\_w\_heat.kfx**
  2. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
5. Once the protocol is completed, the “Lysis Plate” will contain lysate that is ready to proceed to Part 2 (**\*caution\* sample may be hot**).

#### 2. Nanotrap Microbiome Combined NucleoMag KingFisher Apex Procedure-Part 2

1. *Prepare* “NM Binding” Plate
  1. To a new KingFisher 96 Deep Well Plate, add 450  $\mu\text{L}$  of the cleared lysate (NT lysate) from each well of the lysis plate used in “Part 1 step 5” of the protocol. Keep track of which well contains which sample in this new bead binding plate.

2. Add 475  $\mu\text{L}$  of Binding Buffer MWA2 to each well in which lysate was added.
  3. Vortex the NucleoMag B-beads thoroughly and add 25  $\mu\text{L}$  to each well.
    - a) Note: Binding mix (MWA2 + B-beads) can be pre-mixed before their addition to the plate.
  2. Prepare “1<sup>st</sup> Wash MWA3” Plate
    1. Add 850  $\mu\text{L}$  of Wash Buffer MWA3 to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- “NM Binding” Plate wells.
  3. Prepare “2<sup>nd</sup> Wash MWA3” Plate
    1. Add 850  $\mu\text{L}$  of Wash Buffer MWA3 to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- “NM Binding” Plate wells.
  4. Prepare “3<sup>rd</sup> Wash MWA4” Plate
    1. Add 850  $\mu\text{L}$  of Wash Buffer MWA4 to a new KingFisher 96 Deep Well Plate matching the number and location of the KingFisher 96 Deep Well Plate- “NM Binding” Plate wells.
  5. Prepare “Elution” Plate
    1. Add 100  $\mu\text{L}$  of Rnase-free water to a new KingFisher 96- 200  $\mu\text{L}$  plate matching the number and location of the KingFisher 96 Deep Well Plate- “NM Binding” Plate wells.
  6. Prepare “Tip Plate”
    1. Insert the KingFisher 96 Deep Well Comb into a new KingFisher 96 Deep Well Plate
  7. *Run Extraction Kit Script (Request file at sales@ceresnano.com)*
    1. Run  
**NucleoMag\_DNA\_RNA\_Water\_CeresNanoTrap\_Apex\_Rev02.kfx**
    2. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
3. Once the protocol is completed, the KingFisher 96-Elution Plate contains eluates that are ready for downstream analysis or can be stored at  $-80^{\circ}\text{C}$ .  
*Note: Multiple freeze-thaw cycles may cause degradation.*

---

**Attachments: 2**

*KingFisher™ Apex*

1. *NT\_Microbiome\_A\_and\_B\_NucleoMag®\_24\_w\_heat.kfx*
2. *NucleoMag\_DNA\_RNA\_Water\_CeresNanoTrap\_Apex\_Rev02.kfx*