

Automated Nanotrap[®] Wastewater Protocol using Promega Maxwell[®] HT Environmental TNA Kit

Objective: This protocol uses Nanotrap[®] Magnetic Virus Particles and Nanotrap[®] Enhancement Reagent 2 to capture and concentrate viruses in wastewater samples. It is optimized for viral capture from 10 mL samples of wastewater and is compatible with Promega Maxwell[®] HT Environmental TNA Kit. This method has been verified with SARS-CoV-2 viral samples.

Materials and equipment:

1. Wastewater sample
2. Nanotrap[®] Enhancement Reagent 2 (ER2) (Ceres Nanosciences SKU 10112)
 - a. Using ER2 improves viral detection by 1-2 Ct values when used with Nanotrap[®] Magnetic Virus Particles, but it is ok to skip addition of ER2 in this protocol if it is not available.
3. Nanotrap[®] Magnetic Virus Particles (Ceres Nanosciences SKU 44202)
4. Thermo Scientific KingFisher[™] Apex
5. Promega Maxwell[®] HT Environmental TNA Kit (Promega Cat# AX9190)
 - a. Wash Buffer
 - b. Resin
 - c. Tris-HCl (pH 8.0)
6. 5.5M Guanidine Thiocyanate (Sigma Cat# G9277-100G)
7. Isopropanol
8. 50% Ethanol
9. 3-24 well-KF Deep Well Plates
10. 1-24 well-KF Deep Well Comb
11. 4-96 well-KF Deep Well Plates
12. 1-96 well-KF 200 μ L Micro Well Plate
13. 1-96 well-KF Deep Well Comb

Procedure:

Capture and Concentrate Virus using Nanotrap[®] Particles

1. Ceres Nanotrap KF Procedure - Part 1

a. Prepare Sample Plates

- i. Add 4,875 μ L of wastewater sample from wastewater bottle to one well (one well per sample) of a new KingFisher[™] 24 Well Deep Well Plate.
- ii. Add another 4,875 μ L of each sample to the same well on a second KingFisher[™] 24 Well Deep Well Plate. For example, if you loaded a sample into well A1 of the first plate, load the second volume of that sample into well A1 of the second plate.
- iii. Add 50 μ L of Nanotrap[®] Enhancement Reagent 2 (ER2) Solution to each wastewater sample on the two KingFisher[™] 24 Well Deep Well sample plates.
 1. Add 75 μ L of Nanotrap[®] Magnetic Virus Particles to each wastewater sample on the two KingFisher[™] 24 Well Deep Well sample plates.
- iv. Insert the comb into one of the sample plates. This will be “Sample Plate 1” while the other plate will be “Sample Plate 2”.

Using ER2 improves viral detection by 1-2 Ct values when used with Nanotrap[®] Magnetic Virus Particles, but it is ok to skip addition of ER2 and proceed to the next step.

- b. *Prepare Lysis Plate*
 - i. Add 300 μ L of guanidine thiocyanate solution to the third Kingfisher[™] 24 Well Deep Well Plate matching the number and location of the “Sample Plate” wells.
- c. *Run NT KingFisher[™] Protocol (See attached file)*
 - i. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
- d. Once the protocol is completed, the KingFisher[™] plate to which the lysis solution was added will contain lysate that is ready to be purified on the KingFisher[™] Apex.

2. Promega Maxwell[®] KF Extraction Procedure -Part 2

- a. *Prepare Wash Plate 1*
 - i. Add 900 μ L of Wash buffer + 100 μ L of 50% Ethanol to a new KingFisher[™] 96 Deep Well Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- b. *Prepare Wash Plate 2*
 - i. Add 900 μ L of Wash buffer + 100 μ L of 50% Ethanol to a new KingFisher[™] 96 Deep Well Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- c. *Prepare Wash Plate 3*
 - i. Add 450 μ L of 50% Ethanol to a new KingFisher[™] 96 Deep Well Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- d. *Prepare Elution Plate*
 - i. Add 50 μ L of Tris-HCL Elution buffer to a new KingFisher[™] 96- 200 μ L plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- e. *Prepare Promega Bead Binding Plate*
 - i. To a new KingFisher[™] 96 Deep Well Plate, add 300 μ L of the lysate from each well of the lysis plate used in Part 1 of the protocol. Keep track of which well contains which sample in this new bead binding plate.
 - ii. Add 400 μ L of Isopropanol to each well in which lysate was added.
 - iii. Vortex the Resin Beads thoroughly to homogenize and add 35 μ L to each well.
 - iv. Insert the KingFisher[™] 96 Deep Well Comb into the bead binding plate.
- f. *Run Promega KingFisher[™] Protocol (See attached file)*
 - i. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.

3. Once the protocol is completed, the KingFisher[™] 96 -Elution Plate will contain purified viral RNA that is ready to be loaded onto a PCR plate.

Attachments: 2

KingFisher[™] Apex

Using ER2 improves viral detection by 1-2 Ct values when used with Nanotrap[®] Magnetic Virus Particles, but it is ok to skip addition of ER2 and proceed to the next step.

1. *KF-005-WW-Nanotrap-24.kfx*
2. *KF-005-WW-Promega-96.kfx*