

Automated Nanotrap[®] Wastewater Protocol using ZymoBIOMICS[®] MagBead DNA/RNA Kit

Objective: This protocol uses Nanotrap[®] Magnetic Virus Particles and Nanotrap[®] Enhancement Reagent 2 to capture and concentrate viruses in wastewater samples. It is optimized for viral capture from 10 mL samples of wastewater and is compatible with ZymoBIOMICS[®] MagBead DNA/RNA Kit. This method has been verified with SARS-CoV-2 viral samples.

Materials and equipment:

1. Wastewater sample
2. Nanotrap[®] Enhancement Reagent 2 (ER2) (Ceres Nanosciences SKU 10112)
 - a. Using ER2 improves viral detection by 1-2 Ct values when used with Nanotrap[®] Magnetic Virus Particles, but it is ok to skip addition of ER2 in this protocol if it is not available.
3. Nanotrap[®] Magnetic Virus Particles (Ceres Nanosciences SKU 44202)
4. ZymoBIOMICS[®] MagBead DNA/RNA Kit (Zymo Cat# R2135)
 - a. DNA/RNA Lysis Buffer
 - b. MagBinding Beads
 - c. Wash Buffer 1
 - d. Wash Buffer 2
 - e. DNase/RNase Free Water
5. Thermo Scientific KingFisher[™] Apex
6. Molecular Grade Ethanol
7. 3-24 well-KF Deep Well Plates
8. 1-24 well-KF Deep Well Comb
9. 4-96 well-KF Deep Well Plates
10. 1-96 well-KF Plate (200 μ L)
11. 1-96 well-KF Deep Well Comb

Procedure:

Capture and Concentrate Virus using Nanotrap[®] Particles

1. Ceres Nanotrap KF Procedure - Part 1

a. Prepare Sample Plates

- i. Add 4,875 μ L of wastewater sample from wastewater bottle to one well (one well per sample) of a new KingFisher[™] 24 Well Deep Well Plate.
- ii. Add another 4,875 μ L of each sample to the same well on a second KingFisher[™] 24 Well Deep Well Plate. For example, if you loaded a sample into well A1 of the first plate, load the second volume of that sample into well A1 of the second plate.
- iii. Add 50 μ L of Nanotrap[®] Enhancement Reagent 2 (ER2) to each wastewater sample on the two KingFisher 24 Well Deep Well sample plates.
 1. Using ER2 improves viral detection by 1-2 Ct values when used with Nanotrap[®] Magnetic Virus Particles, but it is ok to skip addition of ER2 and proceed to the next step.
- iv. Add 75 μ L of Nanotrap[®] Magnetic Virus Particles to each wastewater sample on the two KingFisher 24 Well Deep Well sample plates.

- v. Insert the comb into one of the sample plates. This will be “Sample Plate 1” while the other plate will be “Sample Plate 2”.
- b. *Prepare Lysis Plate*
 - i. Add 400 μL of Zymo Lysis Buffer to the third KingFisher™ 24 Well Deep Well Plate matching the number and location of the “Sample Plate” wells.
- c. *Run NT KingFisher™ Protocol (See attached file)*
 - i. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
- d. Once the protocol is completed, the KingFisher™ plate to which a lysis solution was added will contain lysate that is ready to be run on the KingFisher™ Zymo MagBead extraction protocol.

2. ZymoBIOMICS MagBead KF Extraction Procedure -Part 2

- a. *Prepare Wash Plate 1*
 - i. Add 500 μL of Wash Buffer 1 to a new KingFisher™ 96 Deep Well Plate matching the number and location of the KingFisher™ 96 Deep Well Plate- Bead Binding Plate wells.
- b. *Prepare Wash Plate 2*
 - i. Add 500 μL of Wash Buffer 2 to a new KingFisher™ 96 Deep Well Plate matching the number and location of the KingFisher™ 96 Deep Well Plate- Bead Binding Plate wells.
- c. *Prepare Wash Plate 3*
 - i. Add 500 μL of Ethanol to a new KingFisher™ 96 Deep Well Plate matching the number and location of the KingFisher™ 96 Deep Well Plate- Bead Binding Plate wells.
- d. *Prepare Wash Plate 4*
 - i. Add 500 μL of Ethanol to a new KingFisher™ 96 Deep Well Plate matching the number and location of the KingFisher™ 96 Deep Well Plate- Bead Binding Plate wells.
- e. *Prepare Elution Plate*
 - i. Add 50 μL of Zymo DNase/RNase-free water to a new KingFisher™ 96-Well (200 μL) plate matching the number and location of the KingFisher™ 96 Deep Well Plate- Bead Binding Plate wells.
- f. *Prepare Zymo Bead Binding Plate*
 - i. To a new KingFisher™ 96 Deep Well Plate, add 400 μL of the lysate from each well of the Lysis Plate used in Part 1 of the protocol. Keep track of which well contains which sample in this new Bead Binding Plate.
 - ii. Add 400 μL of Ethanol to each well in which lysate was added.
 - iii. Vortex the Zymo MagBinding Beads thoroughly and add 30 μL to each well.
 - iv. Insert the KingFisher™ 96 Deep Well Comb into the bead binding plate.
- g. *Run Zymo MagBead KingFisher™ Protocol (See attached file)*
 - i. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.

3. Once the protocol is completed, the KingFisher™ 96 -Elution Plate will contain purified viral RNA that is ready to be loaded onto a PCR plate.

Attachments: 2

KingFisher™ Apex

1. *KF-007-WW-Nanotrap-24.kfx*
2. *KF-007-WW-Zymo-96.kfx*