

Automated Nanotrap[®] Wastewater Protocol using QIAGEN MagAttract[®] Viral RNA Kit

Objective: This protocol uses Nanotrap[®] Magnetic Virus Particles and Nanotrap[®] Enhancement Reagent 2 to capture and concentrate viruses in wastewater samples. It is optimized for viral capture from 10 mL samples of wastewater and is compatible with QIAGEN MagAttract[®] Viral RNA Kit. This method has been verified with SARS-CoV-2 viral samples.

Materials and equipment:

1. Wastewater sample
2. Nanotrap[®] Enhancement Reagent 2 (ER2) (Ceres Nanosciences SKU 10112)
 - a. Using ER2 improves viral detection by 1-2 Ct values when used with Nanotrap[®] Magnetic Virus Particles, but it is ok to skip addition of ER2 in this protocol if it is not available.
3. Nanotrap[®] Magnetic Virus Particles (Ceres Nanosciences SKU 44202)
4. Thermo Scientific KingFisher[™] Apex
5. MagAttract[®] Viral RNA Kit (QIAGEN Cat# 955538)
 - a. Buffer ACL
 - b. Buffer QSB1
 - c. Buffer MW1
 - d. Buffer AVE
 - e. MagAttract Beads
6. 5.5M Guanidine Thiocyanate (Sigma Cat# G9277-100G)
7. 80% Ethanol
8. 3-24 well-KF Deep Well Plates
9. 1-24 well-KF Deep Well Comb
10. 5-96 well-KF Deep Well Plates
11. 1-96 well-KF 200 μ L Well Plate
12. 1-96 well-KF Deep Well Comb

Procedure:

Capture and Concentrate Virus using Nanotrap[®] Particles

1. Ceres Nanotrap KF Procedure - Part 1

a. Prepare Sample Plates

- i. Add 4,875 μ L of wastewater sample from wastewater bottle to one well (one well per sample) of a new KingFisher[™] 24 Well Deep Well Plate.
- ii. Add another 4,875 μ L of each sample to the same well on a second KingFisher[™] 24 Well Deep Well Plate. For example, if you loaded a sample into well A1 of the first plate, load the second volume of that sample into well A1 of the second plate.
- iii. Add 50 μ L of Nanotrap[®] Enhancement Reagent 2 (ER2) Solution to each wastewater sample on the two KingFisher[™] 24 Well Deep Well sample plates.
 1. Using ER2 improves viral detection by 1-2 Ct values when used with Nanotrap[®] Magnetic Virus Particles, but it is ok to skip addition of ER2 and proceed to the next step.

- iv. Add 75 μL of Nanotrap[®] Magnetic Virus Particles to each wastewater sample on the two KingFisher 24 Well Deep Well sample plates.
 - v. Insert the comb into one of the sample plates. This will be “Sample Plate 1” while the other plate will be “Sample Plate 2”.
- b. *Prepare Lysis Plate*
 - i. Add 300 μL of MagAttract[®] Guanidine Thiocyanate solution to the third Kingfisher[™] 24 Well Deep Well Plate matching the number and location of the “Sample Plate” wells.
 - c. *Run NT KingFisher[™] Protocol (See attached file)*
 - i. Follow the on-screen instructions loading the previously prepared plates at the appropriate time.
 - d. Once the protocol is completed, the KingFisher[™] plate to which lysis solution was added will contain lysate that is ready to be run on the Kingfisher’s MagAttract[®] procedure.

2. MagAttract KF Extraction Procedure -Part 2

- a. *Prepare Wash Plate 1*
 - i. Add 500 μL of Buffer MW1 to a new KingFisher[™] 96 Deep Well Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- b. *Prepare Wash Plate 2*
 - i. Add 500 μL of Buffer MW1 to a new KingFisher[™] 96 Deep Well Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- c. *Prepare Wash Plate 3*
 - i. Add 500 μL of 80% Ethanol to a new KingFisher[™] 96 Deep Well Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- d. *Prepare Wash Plate 4*
 - i. Add 500 μL of 80% Ethanol to a new KingFisher[™] 96 Deep Well Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- e. *Prepare Elution Plate*
 - i. Add 50 μL of Buffer AVE to a new KingFisher[™] 96- Well (200 μL) Plate matching the number and location of the KingFisher[™] 96 Deep Well Plate- Bead Binding Plate wells.
- f. *Prepare MagAttract[®] Bead Binding Plate*
 - i. To a new KingFisher[™] 96 Deep Well Plate, add 300 μL of the lysate from each well of the lysis plate used in Part 1 of the protocol. Keep track of which well contains which sample in this new bead binding plate.
 - ii. Add 170 μL of Buffer ACL and 275 μL of Buffer QSB1 to each well in which lysate was added.
 - iii. Vortex the MagAttract[®] Beads thoroughly and add 25 μL to each well.
 - iv. Insert the KingFisher[™] 96 Deep Well Comb into the Bead Binding Plate.

